

Supplementary file

Supplementary Table: Antibodies used for immunohistochemistry.

Antigen	Clone	Source	Dilution	Preparation
Vimentin	V9 (mm)	Cell marque	x200	Autostainer Ventata
AE1AE3	AE1AE3 (mm)	Dako	x300	Autostainer Ventata
S100	(rp)	Dako	x5000	Autostainer Dako
CD34	QBEND/10	Cell marque	x200	Autostainer Ventata
CD3	LN10 (mm)	Novocastra	x300	Autostainer Ventata
CD20	L26 (mm)	Dako	x600	Autostainer Dako
Desmin	D33 (mm)	Dako	x300	Autostainer Ventata
Smooth muscle actin	1A4 (mm)	Dako	x1000	Autostainer Dako
CD68	PGM1 (mm)	Dako	x300	Autostainer Dako
ALK.D5F3	D5F3 (rm)	Ventana	Ready to use, x1	Autostainer Ventata
CD117	(rp)	Thermo scientific	x400	Autostainer Ventata
PLAP	8A9 (mm)	Thermo scientific	x100	Autostainer Ventata
Beta-hCG	(rp)	Gennova	x300	Autostainer Ventata
Oct4	MRQ-10 (mm)	Cell marque	x100	Autostainer Ventata
Alfa-fetoprotein	(rp)	Cell marque	x400	Autostainer Ventata